



Quality of life since 1825

# Washbasin bundle Nautic 60 cm

All In One

## Product features

- Complete package including washbasin, mixer, bottom valve and water trap
- Sealed bracket holes and overflow channel for use in hospitals
- Equipted with elongated lever

## Article numbers

Art-no.	EAN
GB116055609964R1	7391530087153

## Packaging information

Length: 605 mm | Width: 210 mm | Height: 455 mm  
Weight: 18.1 kg  
Weight: 1 pc

## Assembling & Care

Washbasin mounted directly on wall with bolts. Nuts included. Prior to installation the lines must be flushed up to the faucet. The faucet is connected with the cold water line on the right and the hot water line on the left. Water pressure: max 1000 kPa, min 50kPa. Hot water temperature: max 70°C, min 45°C. Please note that we are not liable for damage resulting from improper installation. Clean the porcelain often to prevent dirt and limescale deposits forming. If your regular detergent does not give the desired result, try citric acid or tartaric acid. Wipe/rub the damp surface dry and rinse thoroughly. Remove limescale deposits with household vinegar. Abbrasive cleaning agents must not be used. Never use abrasive detergents. They damage the glaze over time. And never use strongly alkaline, lye-based blockage removers that can damage the porcelain surface and harm the environment.

## Energy labeling and certifications

