



Quality of life since 1825

Armrest Nautic 305502

SQ/QR, white

Product features

- Ergonomically designed grab bar in dark grey contrast colour
- The robust seat buffers make the seat extra sturdy and reliable
- Fits all toilets in the Nautic series, except wall-mounted 1530/5530

Article numbers

Art-no.
GB88305502

EAN
7391530091563

Assembling & Care

Assemble according to the enclosed assembly instructions Use soapy water or a mild cleaning agent. Never use chlorine-based or abrasive cleaners. Then wipe with a clean cloth. Wipe away limescale and urine stains under the seat ring so they don't become ingrained. If the stain doesn't come off, rub repeatedly with warm soapy water using a soft cloth. Do not spray cleaning agents directly onto the seat, as internal components may be damaged. Use a soft cloth for wiping. Hard seats may develop surface wear under the seat ring if cleaning intervals are too long. These can possibly be refreshed with a light polish using car polish without abrasives on a soft cloth. Make it a habit to leave the lid and seat up after cleaning the toilet bowl so that detergent and water have a chance to dry.

Energy labeling and certifications

