



Quality of life since 1825

Toilet seat – Care 9MG7

integrated grab rail, white

Product features

- Fits wall hung toilet 4G01 & 4G95
- Slip stop for side stability
- Integrated grab rail for opening lid

Article numbers

Art-no.
9MG76101

EAN
7391530094212



Packaging information

Length: 500 mm | Width: 462 mm | Height: 59 mm
Weight: 3.35 kg
Weight: 1 pc

Assembling & Care

See the enclosed assembly instructions. Use soapy water or a mild cleaning agent. Never use chlorine-based or abrasive cleaners. Then wipe with a clean cloth. Wipe away limescale and urine stains under the seat ring so they don't become ingrained. If the stain doesn't come off, rub repeatedly with warm soapy water using a soft cloth. Do not spray cleaning agents directly onto the seat, as internal components may be damaged. Use a soft cloth for wiping. Hard seats may develop surface wear under the seat ring if cleaning intervals are too long. These can possibly be refreshed with a light polish using car polish without abrasives on a soft cloth. Make it a habit to leave the lid and seat up after cleaning the toilet bowl so that detergent and water have a chance to dry.

Energy labeling and certifications

