



Quality of life since 1825

Toilet seat Nautic – Care

for armrest, SQ/QR, grey

Product features



- Seat ring that fits both young and adult users
- The robust seat buffers make the seat extra sturdy and reliable
- Soft Close & Quick release seat

Article numbers

Art-no.
8M93S116

EAN
7391530096094

Packaging information

Length: 460 mm | Width: 383 mm | Height: 49 mm

Weight: 2.7 kg

Weight: 1 pc

Assembling & Care

See the enclosed assembly instructions. Use soapy water or a mild cleaning agent. Never use chlorine-based or abrasive cleaners. Then wipe with a clean cloth. Wipe away limescale and urine stains under the seat ring so they don't become ingrained. If the stain doesn't come off, rub repeatedly with warm soapy water using a soft cloth. Do not spray cleaning agents directly onto the seat, as internal components may be damaged. Use a soft cloth for wiping. Hard seats may develop surface wear under the seat ring if cleaning intervals are too long. These can possibly be refreshed with a light polish using car polish without abrasives on a soft cloth. Make it a habit to leave the lid and seat up after cleaning the toilet bowl so that detergent and water have a chance to dry.

Energy labeling and certifications

